




KCNH1 Antibody

Product Code	CSB-PA012035DSR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95259
Immunogen	Recombinant Human Potassium voltage-gated channel subfamily H member 1 protein (840-989AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Potassium voltage-gated channel subfamily H member 1 (Ether-a-go-go potassium channel 1) (EAG channel 1) (h-eag) (hEAG1) (Voltage-gated potassium channel subunit Kv10.1), KCNH1, EAG EAG1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	KCNH1

Image



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA012035DSR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.