





KCNIP2 Antibody

Product Code	CSB-PA889100DSR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NS61
Immunogen	Recombinant Human Kv channel-interacting protein 2 protein (1-90AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kv channel-interacting protein 2 (KChIP2) (A-type potassium channel modulatory protein 2) (Cardiac voltage-gated potassium channel modulatory subunit) (Potassium channel-interacting protein 2), KCNIP2, KCHIP2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	KCNIP2

Image



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA889100DSR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human heart tissue using CSB-PA889100DSR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.