





MYO1C Antibody

Product Code	CSB-PA015340ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00159
Immunogen	Recombinant Human Unconventional myosin-Ic protein (804-1063AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Unconventional myosin-Ic (Myosin I beta) (MMI-beta) (MMIb), MYO1C
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MYO1C

Image



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA015340ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA015340ESR2HU at dilution of 1:100