





RGS7 Antibody

Product Code	CSB-PA019659ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49802
Immunogen	Recombinant Human Regulator of G-protein signaling 7 protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Regulator of G-protein signaling 7 (RGS7), RGS7
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RGS7

Image



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA019659ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human endometrial cancer using CSB-PA019659ESR1HU at dilution of 1:100