




STX5 Antibody

Product Code	CSB-PA619754ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13190
Immunogen	Recombinant Human Syntaxin-5 protein (55-355AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Syntaxin-5, STX5, STX5A
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	STX5

Image



Immunohistochemistry of paraffin-embedded human lung tissue using CSB-PA619754ESR1HU at dilution of 1:100

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.