




TIMM8A Antibody

Product Code	CSB-PA023557ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60220
Immunogen	Recombinant Human Mitochondrial import inner membrane translocase subunit Tim8 A protein (1-97AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitochondrial import inner membrane translocase subunit Tim8 A (Deafness dystonia protein 1) (X-linked deafness dystonia protein), TIMM8A, DDP DDP1 TIM8A
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TIMM8A

Image



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA023557ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.