





TRAF4 Antibody

Product Code	CSB-PA024152ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BUZ4
Immunogen	Recombinant Human TNF receptor-associated factor 4 protein (251-470AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	TNF receptor-associated factor 4 (Cysteine-rich domain associated with RING and Traf domains protein 1) (Metastatic lymph node gene 62 protein) (MLN 62) (RING finger protein 83), TRAF4, CART1 MLN62 RNF83
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TRAF4

Image



Immunohistochemistry of paraffin-embedded human lung cancer using CSB-PA024152ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA024152ESR1HU at dilution of 1:100

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.