




TRIM32 Antibody

Product Code	CSB-PA614389ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13049
Immunogen	Recombinant Human E3 ubiquitin-protein ligase TRIM32 protein (1-280AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase TRIM32 (EC 2.3.2.27) (72 kDa Tat-interacting protein) (RING-type E3 ubiquitin transferase TRIM32) (Tripartite motif-containing protein 32) (Zinc finger protein HT2A), TRIM32, HT2A
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TRIM32

Image



Immunohistochemistry of paraffin-embedded human prostate cancer using CSB-PA614389ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.