





VPS35 Antibody

Product Code	CSB-PA839401ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96QK1
Immunogen	Recombinant Human Vacuolar protein sorting-associated protein 35 protein (1-250AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vacuolar protein sorting-associated protein 35 (hVPS35) (Maternal-embryonic 3) (Vesicle protein sorting 35), VPS35, MEM3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	VPS35

Image



Immunohistochemistry of paraffin-embedded human ovarian cancer using CSB-PA839401ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA839401ESR2HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.