




NME4 Antibody

Product Code	CSB-PA015889ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00746
Immunogen	Recombinant Human Nucleoside diphosphate kinase, mitochondrial protein (1-187AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nucleoside diphosphate kinase, mitochondrial (NDK) (NDP kinase, mitochondrial) (EC 2.7.4.6) (Nucleoside diphosphate kinase D) (NDPKD) (nm23-H4), NME4, NM23D
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NME4

Image



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA015889ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.