





NPFF Antibody

Product Code	CSB-PA015981ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15130
Immunogen	Recombinant Human Pro-FMRFamide-related neuropeptide FF protein (34-113AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pro-FMRFamide-related neuropeptide FF (FMRFamide-related peptides) [Cleaved into: Neuropeptide SF (NPSF); Neuropeptide FF (NPFF); Neuropeptide AF (NPAF)], NPFF
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NPFF

Image



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA015981ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA015981ESR2HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.