





CAPNS1 Antibody

Product Code	CSB-PA004503ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04632
Immunogen	Recombinant Human Calpain small subunit 1 protein (69-268AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calpain small subunit 1 (CSS1) (Calcium-activated neutral proteinase small subunit) (CANP small subunit) (Calcium-dependent protease small subunit) (CDPS) (Calcium-dependent protease small subunit 1) (Calpain regulatory subunit), CAPNS1, CAPN4 CAPNS
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CAPNS1

Image



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA004503ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human gastric cancer using CSB-PA004503ESR2HU at dilution of 1:100

