





CCNE2 Antibody

Product Code	CSB-PA004818ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O96020
Immunogen	Recombinant Human G1/S-specific cyclin-E2 protein (255-404AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	G1/S-specific cyclin-E2, CCNE2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CCNE2

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA004818ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human thyroid tissue using CSB-PA004818ESR2HU at dilution of 1:100