





MMP19 Antibody

Product Code	CSB-PA857450ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99542
Immunogen	Recombinant Human Matrix metalloproteinase-19 protein (209-508AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Matrix metalloproteinase-19 (MMP-19) (EC 3.4.24.-) (Matrix metalloproteinase RASI) (Matrix metalloproteinase-18) (MMP-18), MMP19, MMP18 RASI
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	MMP19

Image



Western blot
 All lanes: MMP19 antibody at 2.58µg/ml + Mouse lung tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 58, 26, 7, 34 kDa
 Observed band size: 58 kDa



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA857450ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.