




CD200 Antibody

Product Code	CSB-PA004895ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41217
Immunogen	Recombinant Human OX-2 membrane glycoprotein protein (30-240AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	OX-2 membrane glycoprotein (CD antigen CD200), CD200, MOX1 MOX2
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CD200

Image



Western blot
 All lanes: CD200antibody at 0.71μg/ml + Mouse kidney tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 32, 31, 33 kDa
 Observed band size: 32 kDa

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.