





CYP4B1 Antibody

Product Code	CSB-PA006448ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13584
Immunogen	Recombinant Human Cytochrome P450 4B1 protein (30-250AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytochrome P450 4B1 (EC 1.14.14.1) (CYP4B1) (Cytochrome P450-HP), CYP4B1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	CYP4B1

Image



Immunohistochemistry of paraffin-embedded human lung cancer using CSB-PA006448ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA006448ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.