




RALB Antibody

Product Code	CSB-PA019297ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11234
Immunogen	Recombinant Human Ras-related protein Ral-B protein (1-206AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ras-related protein Ral-B, RALB
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RALB

Image



Western blot

All lanes: RALB antibody at 1.26µg/ml

Lane 1: Mouse brain tissue

Lane 2: Mouse lung tissue

Lane 3: MCF-7 whole cell lysate

Lane 4: HepG2 whole cell lysate

Lane 5: Jurkat whole cell lysate

Lane 6: NIH/3T3 whole cell lysate


Lane 7: Hela whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 24, 26 kDa

Observed band size: 26 kDa



Immunohistochemistry of paraffin-embedded

human ovarian cancer using CSB-

PA019297ESR2HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.