




# DCTD Antibody

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | CSB-PA006563ESR1HU   |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.  |
| <b>Uniprot No.</b>         | P32321   |
| <b>Immunogen</b>           | Recombinant Human Deoxycytidylate deaminase protein (1-178AA)  |
| <b>Raised In</b>           | Rabbit   |
| <b>Species Reactivity</b>  | Human  |
| <b>Tested Applications</b> | ELISA, IHC; Recommended dilution: IHC:1:20-1:200               |
| <b>Form</b>                | Liquid   |
| <b>Conjugate</b>           | Non-conjugated   |
| <b>Storage Buffer</b>      | PBS with 0.02% sodium azide, 50% glycerol, pH7.3.              |
| <b>Purification Method</b> | Antigen Affinity Purified                                      |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Alias</b>               | Deoxycytidylate deaminase (EC 3.5.4.12) (dCMP deaminase), DCTD |
| <b>Immunogen Species</b>   | Homo sapiens (Human)   |
| <b>Research Area</b>       | Epigenetics and Nuclear Signaling                              |
| <b>Target Names</b>        | DCTD   |

## Image



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA006563ESR1HU at dilution of 1:100

**Usage** For Research Use Only. Not for use in diagnostic or therapeutic procedures.