




DCTD Antibody

Product Code	CSB-PA006563ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P32321
Immunogen	Recombinant Human Deoxycytidylate deaminase protein (1-178AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Deoxycytidylate deaminase (EC 3.5.4.12) (dCMP deaminase), DCTD
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	DCTD

Image


Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA006563ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.