





RASGRF1 Antibody

Product Code	CSB-PA618811ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13972
Immunogen	Recombinant Human Ras-specific guanine nucleotide-releasing factor 1 protein (1-190AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ras-specific guanine nucleotide-releasing factor 1 (Ras-GRF1) (Guanine nucleotide-releasing protein) (GNRP) (Ras-specific nucleotide exchange factor CDC25), RASGRF1, CDC25 GNRP GRF1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RASGRF1

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA618811ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA618811ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.