





RPLP2 Antibody

Product Code	CSB-PA020342DSR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05387
Immunogen	Recombinant Human 60S acidic ribosomal protein P2 protein (1-115AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	60S acidic ribosomal protein P2 (Large ribosomal subunit protein P2) (Renal carcinoma antigen NY-REN-44), RPLP2, D11S2243E RPP2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	RPLP2

Image



Immunohistochemistry of paraffin-embedded human lung tissue using CSB-PA020342DSR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA020342DSR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.