





ABCA3 Antibody

Product Code	CSB-PA859942ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99758
Immunogen	Recombinant Human ATP-binding cassette sub-family A member 3 protein (43-203AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	ATP-binding cassette sub-family A member 3 (ABC-C transporter) (ATP-binding cassette transporter 3) (ATP-binding cassette 3), ABCA3, ABC3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ABCA3

Image



Immunohistochemistry of paraffin-embedded human lung tissue using CSB-PA859942ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA859942ESR1HU at dilution of 1:100