




AP3B1 Antibody

Product Code	CSB-PA001874ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00203
Immunogen	Recombinant Human AP-3 complex subunit beta-1 protein (895-1094AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	AP-3 complex subunit beta-1 (Adaptor protein complex AP-3 subunit beta-1) (Adaptor-related protein complex 3 subunit beta-1) (Beta-3A-adaptin) (Clathrin assembly protein complex 3 beta-1 large chain), AP3B1, ADTB3A
Immunogen Species	Homo sapiens (Human)
Research Area	Microbiology
Target Names	AP3B1

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA001874ESR1HU at dilution of 1:100