





ATP2B2 Antibody

Product Code	CSB-PA002336DSR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01814
Immunogen	Recombinant Human Plasma membrane calcium-transporting ATPase 2 protein (1-80AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Plasma membrane calcium-transporting ATPase 2 (PMCA2) (EC 3.6.3.8) (Plasma membrane calcium ATPase isoform 2) (Plasma membrane calcium pump isoform 2), ATP2B2, PMCA2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ATP2B2

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA002336DSR1HU at dilution of 1:100



Immunofluorescent analysis of 293 cells using CSB-PA002336DSR1HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.