





HLA-DOA Antibody

Product Code	CSB-PA356784ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06340
Immunogen	Recombinant Human HLA class II histocompatibility antigen, DO alpha chain protein (26-135AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	HLA class II histocompatibility antigen, DO alpha chain (MHC DN-alpha) (MHC DZ alpha) (MHC class II antigen DOA), HLA-DOA, HLA-DNA HLA-DZA
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	HLA-DOA

Image



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using CSB-PA356784ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human skin tissue using CSB-PA356784ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.