





NFYB Antibody

Product Code	CSB-PA015774ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25208
Immunogen	Recombinant Human Nuclear transcription factor Y subunit beta protein (1-207AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nuclear transcription factor Y subunit beta (CAAT box DNA-binding protein subunit B) (Nuclear transcription factor Y subunit B) (NF-YB), NFYB, HAP3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NFYB

Image



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA015774ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human thymus tissue using CSB-PA015774ESR2HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.