





PDCD5 Antibody

Product Code	CSB-PA017671ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14737
Immunogen	Recombinant Human Programmed cell death protein 5 protein (1-125AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Programmed cell death protein 5 (TF-1 cell apoptosis-related protein 19) (Protein TFAR19), PDCD5, TFAR19
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PDCD5

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA017671ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA017671ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.