





PDHB Antibody

Product Code	CSB-PA017717ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11177
Immunogen	Recombinant Human Pyruvate dehydrogenase E1 component subunit beta, mitochondrial protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pyruvate dehydrogenase E1 component subunit beta, mitochondrial (PDHE1-B) (EC 1.2.4.1), PDHB, PHE1B
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	PDHB

Image



Immunohistochemistry of paraffin-embedded human gastric cancer using CSB-PA017717ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human heart tissue using CSB-PA017717ESR2HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.