




# RASSF5 Antibody

<b>Product Code</b>	CSB-PA819906ESR1HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q8WWW0
<b>Immunogen</b>	Recombinant Human Ras association domain-containing protein 5 protein (1-265AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Ras association domain-containing protein 5 (New ras effector 1) (Regulator for cell adhesion and polarization enriched in lymphoid tissues) (RAPL), RASSF5, NORE1 RAPL
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Cell Biology
<b>Target Names</b>	RASSF5

## Image



### Western blot

All lanes: RASSF5antibody at 5.46µg/ml

Lane 1: Thp-1 whole cell lysate

Lane 2: 293T whole cell lysate

Lane 3: Jurkat whole cell lysate


Lane 4: A549 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 48, 31, 38 kDa

Observed band size: 48 kDa



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA819906ESR1HU at dilution of 1:100

---

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.