




RASSF5 Antibody

Product Code	CSB-PA819906ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WWW0
Immunogen	Recombinant Human Ras association domain-containing protein 5 protein (1-265AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ras association domain-containing protein 5 (New ras effector 1) (Regulator for cell adhesion and polarization enriched in lymphoid tissues) (RAPL), RASSF5, NORE1 RAPL
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RASSF5

Image



Western blot

All lanes: RASSF5antibody at 5.46μg/ml

Lane 1: Thp-1 whole cell lysate

Lane 2: 293T whole cell lysate

Lane 3: Jurkat whole cell lysate


Lane 4: A549 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 48, 31, 38 kDa

Observed band size: 48 kDa



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA819906ESR1HU at dilution of 1:100