





RDH5 Antibody

Product Code	CSB-PA856423ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92781
Immunogen	Recombinant Human 11-cis retinol dehydrogenase protein (24-318AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	11-cis retinol dehydrogenase (11-cis RDH) (11-cis RoDH) (EC 1.1.1.315) (9-cis retinol dehydrogenase) (9cRDH) (Retinol dehydrogenase 5) (Short chain dehydrogenase/reductase family 9C member 5), RDH5, HSD17B9 RDH1 SDR9C5
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RDH5

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA856423ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human colon tissue using CSB-PA856423ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.