





BFAR Antibody

Product Code	CSB-PA002678ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NZS9
Immunogen	Recombinant Human Bifunctional apoptosis regulator protein (1-140AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Bifunctional apoptosis regulator (RING finger protein 47), BFAR, BAR RNF47
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	BFAR

Image



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA002678ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA002678ESR1HU at dilution of 1:100