





ELMOD1 Antibody

Product Code	CSB-PA839788ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N336
Immunogen	Recombinant Human ELMO domain-containing protein 1 protein (185-334AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	ELMO domain-containing protein 1, ELMOD1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ELMOD1

Image



Immunohistochemistry of paraffin-embedded human lung tissue using CSB-PA839788ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA839788ESR2HU at dilution of 1:100