





EVA1A Antibody

Product Code	CSB-PA862042ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H8M9
Immunogen	Recombinant Human Protein eva-1 homolog A protein (58-152AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein eva-1 homolog A (Protein FAM176A) (Transmembrane protein 166), EVA1A, FAM176A TMEM166
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	EVA1A

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA862042ESR2HU at dilution of 1:100



Immunofluorescent analysis of MCF-7 cells using CSB-PA862042ESR2HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)