





BEX3 Antibody

Product Code	CSB-PA015781ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q00994
Immunogen	Recombinant Human Protein BEX3 protein (1-111AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein BEX3 (Brain-expressed X-linked protein 3) (Nerve growth factor receptor-associated protein 1) (Ovarian granulosa cell 13.0 kDa protein HGR74) (p75NTR-associated cell death executor), BEX3, DXS6984E NADE NGFRAP1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	BEX3

Image



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA015781ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human gastric cancer using CSB-PA015781ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.