




PECR Antibody

Product Code	CSB-PA017769ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BY49
Immunogen	Recombinant Human Peroxisomal trans-2-enoyl-CoA reductase protein (24-303AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peroxisomal trans-2-enoyl-CoA reductase (TERP) (EC 1.3.1.38) (2,4-dienoyl-CoA reductase-related protein) (DCR-RP) (HPDHase) (Short chain dehydrogenase/reductase family 29C member 1) (pVI-ARL), PECR, SDR29C1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	PECR

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA017769ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.