




TNFAIP8 Antibody

Product Code	CSB-PA023960ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95379
Immunogen	Recombinant Human Tumor necrosis factor alpha-induced protein 8 protein (1-198AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tumor necrosis factor alpha-induced protein 8 (TNF alpha-induced protein 8) (Head and neck tumor and metastasis-related protein) (MDC-3.13) (NF-kappa-B-inducible DED-containing protein) (NDED) (SCC-S2) (TNF-induced protein GG2-1), TNFAIP8
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	TNFAIP8

Image



Western blot

All lanes: TNFAIP8 antibody at 2.79μg/ml

Lane 1: A431 whole cell lysate

Lane 2: PC-3 whole cell lysate

Lane 3: Human placenta tissue


Lane 4: Mouse thymus tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 24, 23, 22, 25 kDa

Observed band size: 22 kDa



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA023960ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.