





TRIAP1 Antibody

Product Code	CSB-PA024443ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43715
Immunogen	Recombinant Human TP53-regulated inhibitor of apoptosis 1 protein (1-76AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	TP53-regulated inhibitor of apoptosis 1 (Protein 15E1.1) (WF-1) (p53-inducible cell-survival factor) (p53CSV), TRIAP1, 15E1.1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TRIAP1

Image



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA024443ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA024443ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.