




TRIM69 Antibody

Product Code	CSB-PA774824ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86WT6
Immunogen	Recombinant Human E3 ubiquitin-protein ligase TRIM69 protein (1-341AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase TRIM69 (EC 2.3.2.27) (RFP-like domain-containing protein trimless) (RING finger protein 36) (RING-type E3 ubiquitin transferase TRIM69) (Tripartite motif-containing protein 69), TRIM69, RNF36
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TRIM69

Image



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA774824ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.