




DPF1 Antibody

Product Code	CSB-PA835697ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92782
Immunogen	Recombinant Human Zinc finger protein neuro-d4 protein (261-380AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein neuro-d4 (BRG1-associated factor 45B) (BAF45B) (D4, zinc and double PHD fingers family 1), DPF1, BAF45B NEUD4
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	DPF1

Image



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA835697ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.