




ELF5 Antibody

Product Code	CSB-PA891540ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UKW6
Immunogen	Recombinant Human ETS-related transcription factor Elf-5 protein (1-170AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	ETS-related transcription factor Elf-5 (E74-like factor 5) (Epithelium-restricted ESE-1-related Ets factor) (Epithelium-specific Ets transcription factor 2) (ESE-2), ELF5, ESE2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ELF5

Image



Immunohistochemistry of paraffin-embedded human prostate cancer using CSB-PA891540ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.