




IKBIP Antibody

Product Code	CSB-PA383678ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q70UQ0
Immunogen	Recombinant Human Inhibitor of nuclear factor kappa-B kinase-interacting protein (121-350AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Inhibitor of nuclear factor kappa-B kinase-interacting protein (I kappa-B kinase-interacting protein) (IKBKB-interacting protein) (IKK-interacting protein), IKBIP, IKIP
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	IKBIP

Image



Western blot

All lanes: IKBIPantibody at 2.23μg/ml

Lane 1: Mouse kidney tissue


Lane 2: HL60 whole cell lysate

Secondary


Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 40, 28, 8, 44 kDa

Observed band size: 40, 36 kDa



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA383678ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA383678ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.