





SHOX2 Antibody

Product Code	CSB-PA021277ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60902
Immunogen	Recombinant Human Short stature homeobox protein 2 protein (72-331AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Short stature homeobox protein 2 (Homeobox protein Og12X) (Paired-related homeobox protein SHOT), SHOX2, OG12X SHOT
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	SHOX2

Image



Immunohistochemistry of paraffin-embedded human heart tissue using CSB-PA021277ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human lung tissue using CSB-PA021277ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.