




ADPGK Antibody

Product Code	CSB-PA001380ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BRR6
Immunogen	Recombinant Human ADP-dependent glucokinase protein (218-497AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	ADP-dependent glucokinase (ADP-GK) (ADPGK) (EC 2.7.1.147) (RbBP-35), ADPGK
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ADPGK

Image



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA001380ESR1HU at dilution of 1:100