




ARIH2 Antibody

Product Code	CSB-PA002069ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95376
Immunogen	Recombinant Human E3 ubiquitin-protein ligase ARIH2 protein (194-493AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase ARIH2 (ARI-2) (Protein ariadne-2 homolog) (EC 2.3.2.27) (RING-type E3 ubiquitin transferase ARIH2) (Triad1 protein), ARIH2, ARI2 TRIAD1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ARIH2

Image



Immunohistochemistry of paraffin-embedded human lung tissue using CSB-PA002069ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.