





DNAJB9 Antibody

Product Code	CSB-PA866205ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBS3
Immunogen	Recombinant Human DnaJ homolog subfamily B member 9 protein (24-223AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	DnaJ homolog subfamily B member 9 (Endoplasmic reticulum DNA J domain-containing protein 4) (ER-resident protein ERdj4) (ERdj4) (Microvascular endothelial differentiation gene 1 protein) (Mdg-1), DNAJB9, MDG1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DNAJB9

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA866205ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA866205ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.