




FLVCR2 Antibody

Product Code	CSB-PA892364ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UPI3
Immunogen	Recombinant Human Feline leukemia virus subgroup C receptor-related protein 2 protein (1-100AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Feline leukemia virus subgroup C receptor-related protein 2 (Calcium-chelate transporter) (CCT), FLVCR2, C14orf58
Immunogen Species	Homo sapiens (Human)
Research Area	Microbiology
Target Names	FLVCR2

Image



Western blot

All lanes: FLVCR2 antibody at 0.89μg/ml

Lane 1: MCF-7 whole cell lysate

Lane 2: A549 whole cell lysate

Lane 3: Mouse lung tissue


Lane 4: Mouse liver tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 58, 36 kDa

Observed band size: 58, 86, 45, 26 kDa



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA892364ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.