





MEOX1 Antibody

Product Code	CSB-PA013695ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50221
Immunogen	Recombinant Human Homeobox protein MOX-1 protein (1-254AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Homeobox protein MOX-1 antibody; MEOX1 antibody; MEOX1_HUMAN antibody; Mesenchyme homeo box 1 antibody; Mesenchyme homeobox 1 antibody; MOX 1 antibody; MOX1 antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MEOX1

Image



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA013695ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human gastric cancer using CSB-PA013695ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.