





NLE1 Antibody

Product Code	CSB-PA889125ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NVX2
Immunogen	Recombinant Human Notchless protein homolog 1 protein (206-485AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Notchless protein homolog 1, NLE1
Immunogen Species	Homo sapiens (Human)
Research Area	Stem Cells
Target Names	NLE1

Image



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA889125ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human gastric cancer using CSB-PA889125ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.