




PEX14 Antibody

Product Code	CSB-PA017800ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75381
Immunogen	Recombinant Human Peroxisomal membrane protein PEX14 protein (118-377AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peroxisomal membrane protein PEX14 (PTS1 receptor-docking protein) (Peroxin-14) (Peroxisomal membrane anchor protein PEX14), PEX14
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	PEX14

Image



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA017800ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.